

Intended Use


Thank you for the purchase of our sound amplifier. If you find it's difficult to hearing very well from others a sound amplifier can help you regain your hearing so you can live the life you're used to.

Wearing a sound amplifier means rediscovering the pleasure of gathering with friends and family feeling more confident at work, and enjoying movies, music and TV programs.


Specification

1. Max sound output: 123+4dB
2. Sound gain: ≥ 30 dB
3. Total harmonic wave distortion: $\leq 10\%$
4. Frequency range: 300Hz~4000Hz
5. Input noise: ≤ 32 dB
6. Voltage: D.C. 1.5V
7. Current: ≤ 4 mA

Instruction Diagram



Battery: AG5 or A13 Zinc Air Battery

- "N" is power on
- "OFF" is power off
- "+" Positive side of Battery
- "-" Negative side of Battery
- " " is Volume grade

Operation

1. Open the battery compartment, put into the corresponding battery and close the battery compartment tightly.
2. Insert the earplug to the sound amplifier, put the sound amplifier in your ear and move it up and down to insure it on the right position.
3. Turn on the sound amplifier, then adjust the volume to get clear sound.
4. If you feel the whistling sound, please make sure if you put the earplug tightly, or if you use an appropriate size earplug.



Operation

5. Please take out the battery or turn off the sound amplifier when you stop to use the sound amplifier.
(Note: Be sure the sound amplifier is in "o" position or volume is in the minimum volume of "1", otherwise the instant sound will make you feel uncomfortable)

Maintenance

The sound amplifier doesn't work.

Check as follows:

1. Sound amplifier is power off
2. The battery is in correct direction
3. Run out of battery
4. The earphone is damaged
5. The earphone is broken

Sound is very light:

1. Volume is too low, please increase the volume
2. Battery is low, please replace a new battery

Clean off dirt of the sound amplifier case using a soft dry cloth

Note

1. Please don't use the sound amplifier at a high volume for a long time. It will damage your hearing.
2. Please maintain the product all the time in order to avoid the earplug and microphone opening being blocked by foreign objects.
3. Please don't use when showering, swimming, raining or wet environment, please don't expose the sound amplifier under sunshine, get close to a stove or other heat sources.
4. Please protect the sound amplifier carefully and avoid dropping it on the floor.

Note

5. Please don't remove old batteries in the sound amplifier for a long time in order to avoid possible erosion caused by battery leakage.
6. Please handle the old batteries environmentally friendly.
7. Keep the sound amplifier away from the children.

Sound Amplifier

Operation Manual